Spis treści

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The following Tables are provided in accordance with IEC/ EN 60601-1-2:

Guidance and manufacturer's declaration – electromagnetic emissions		
The Aeroneb Go nebulizer system is intended for use in the electromagnetic environment specified below. The customer or the user of the Aeroneb Go nebulizer system should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The Aeroneb Go nebulizer system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Aeroneb Go nebulizer system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Recommended separation distances between portable and mobile RF communication equipment and the Aeroneb Go

The Aeroneb Go nebulizer system is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Aeroneb Go can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Aeroneb Go nebulizer system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = [1.17] \sqrt{P}$	$d = [1.17] \sqrt{P}$	$d = [2.33] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.75
1	1.17	1.17	2.33
10	3.70	3.70	7.36
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The Aeroneb Go nebulizer system is intended for use in the electromagnetic environment specified below. The customer or the user of the Aeroneb Go should assure that it is used in such an environment.

Immunity test	IEC/EN 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % Ut (>95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles <5 % Ut (>95 % dip in Ut) for 5 sec	<5 % Ut (>95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles <5 % Ut (>95 % dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Aeroneb Go requires continued operation during power mains interruption, it is recommended that the Aeroneb Go be powered from an uninterruptible power supply or battery.

Immunity test	IEC/EN 60601 Test level	Compliance level	Electromagnetic environment – guidance
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: Ut is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The Aeronex Go nebulizer system is intended for use in the electromagnetic environment specified below. The customer or the user of the Aeronex Go should assure that it is used in such an environment.			
Immunity test	IEC/EN 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	[3]V	Portable and mobile RF communications equipment should be used no closer to any part of the Aeronex Go, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [1.17] \sqrt{P}$ $d = [1.17] \sqrt{P} \dots 80\text{MHz to } 800\text{MHz}$ $d = [2.33] \sqrt{P} \dots 800\text{MHz to } 2.5\text{GHz}$
Radiated RF IEC 61000-4-3	3 Vrms 80 MHz to 2.5 GHz	[3]V	where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).

Guidance and manufacturer's declaration – electromagnetic immunity

The Aeroneb Go nebulizer system is intended for use in the electromagnetic environment specified below. The customer or the user of the Aeroneb Go should assure that it is used in such an environment.

Immunity test	IEC/EN 60601 Test level	Compliance level	Electromagnetic environment – guidance
			<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Aeroneb Go is used exceeds the applicable RF compliance level above, the Aeroneb Go should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the Aeroneb Go nebulizer system.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [3]V/m

Diese Seite wurde freigelassen
Cette page est laissée vierge intentionnellement
Deze pagina is met opzet blanco gelaten
Questa pagina è stata lasciata intenzionalmente
in bianco
Ta strona została celowo pozostawiona pusta

Diese Seite wurde freigelassen
Cette page est laissée vierge intentionnellement
Deze pagina is met opzet blanco gelaten
Questa pagina è stata lasciata intenzionalmente
in bianco
Ta strona została celowo pozostawiona pusta

Aerogen Limited,
Galway Business Park,
Dangan,
Galway,
Ireland.
www.aerogen.com
Ph: +353 91 540400

©2013 Aerogen Limited
AG-AG4120-NE Rev C 2013
Manufacturing no.: 30-554 Rev C